

Ligand Optimisation using R

The purpose of these notes is to explain how to apply McGuigan's new extended Ligand Optimisation Method (LOME) using software written in the package R.

The R package

The user is required to download and install the free package R from <http://cran.r-project.org>, if necessary.

Once the user has opened R, running the code involves copying and pasting some commands in a script file into the R window and pressing ENTER, as explained below. It is not necessary for the user to know anything about R itself.


Also, download and install RStudio from:

<https://www.rstudio.com/products/rstudio/download/>

We will use RStudio in the following demonstration. The file.choose command used below works better in RStudio in that the working directory is made available to the user.

The Data Files

Let's assume that the folder CaLOME has been placed on the user's Desktop, as shown below.



The folder CaLOME contains the files that are necessary to run the program LOME on the selected input data. There are two sets of files - one set for the Aequorin data and one set for the fluorescence data - with the contents shown in screenshots.

Aequorin Data Files

calA.csv	the data obtained using the calibrating solutions
bufA.csv	the data obtained from the buffer solutions
parA.csv	the parameter values for nominal ligand concentration and the type of data

The screenshot shows three CSV files open in a text editor. The first two, calA.csv and bufA.csv, contain 14 rows of data. The third, parA.csv, contains 4 rows of parameter values.

File	Row	Content
calA.csv	1	"XT","pX"
	2	2,2.699
	3	1.5,2.8239
	4	1,3
	5	0.5,3.301
	6	0.25,3.6021
	7	0.1,4
	8	0.05,4.301
	9	0.025,4.6021
	10	0.01,5
	11	0.005,5.301
	12	0.0025,5.6021
	13	0.001,6
	14	
bufA.csv	1	"XT","pX"
	2	0.875,5.04435072709255
	3	0.8571,5.2104263381279
	4	0.8333,5.38151259860365
	5	0.8,5.56019378860667
	6	0.75,5.75708605318107
	7	0.6666,5.99658255479781
	8	0.5,6.35398859040876
	9	0.3333,6.68089394069022
	10	0.2,6.99439055463972
	11	0.1,7.35359627377693
	12	
	13	
	14	
parA.csv	1	"Name","Value"
	2	"NomLigConc",1
	3	"Type",2
	4	

Fluorescence Data Files

calF.csv	the data obtained using the calibrating solutions
bufF.csv	the data obtained from the buffer solutions
parF.csv	the parameter values for nominal ligand concentration and the type of data


The screenshot shows three CSV files open in a text editor. The first two, calF.csv and bufF.csv, contain 14 rows of data. The third, parF.csv, contains 4 rows of parameter values.

File	Row	Content
calF.csv	1	"XT","pX"
	2	2,2.699
	3	1.5,2.8239
	4	1,3
	5	0.5,3.301
	6	0.25,3.6021
	7	0.1,4
	8	0.05,4.301
	9	0.025,4.6021
	10	0.01,5
	11	0.005,5.301
	12	0.0025,5.6021
	13	0.001,6
	14	
bufF.csv	1	"XT","pX"
	2	3.5,5.7597
	3	3.4286,5.8605
	4	3.3333,5.9712
	5	3.2,6.0976
	6	3.6,2503
	7	2.6667,6.4525
	8	2.6,7782
	9	1.3333,7.0911
	10	0.8,7.398
	11	0.4,7.7533
	12	
	13	
	14	
parF.csv	1	"Name","Value"
	2	"NomLigConc",4
	3	"Type",1
	4	

These files are in comma separated value (csv) format. A simple way to make data files in this format is to use Excel and then save the file in csv format. Alternatively, a text editor could be used to place the data in csv format; in each row the entries must be separated by commas rather than white space.

Running the LOME program


1. Launch the application RStudio. This will open a window into which R commands can be entered.
2. Open the file LOMEscript.txt. This contains the the required R commands. The lines beginning with `#` are comment lines are are ignored by the R interpreter. Now, both windows will be open on the desktop, as shown below.



3. Set the working directory to tell RStudio where the data files are; output will also be sent to this folder. For this illustration it has been assumed that the files are in a folder named CaLOME and it is placed on the Desktop. The user can place the files in a folder with a name of their choosing and place the folder anywhere in their Documents folder; if so, the `setwd()` command needs editing.

NB. The form of the folder name given here is suitable for a Mac, but in Windows the slashes would be different.


Select `setwd("~/Desktop/CaLOME")` and copy this to the R window at the command prompt `>`, and then press ENTER to run this command.



4. Select the line `source("LOME.txt")` and copy and paste this into the RStudio window at the command prompt, and then press ENTER to run this command. This loads the program LOME into RStudio and this program will be run below once we have entered the data files. In what follows we will enter the data files for an application with Aequorin data.
5. Enter the data files one by one. Select the line

`cacal1 <- read.csv(file.choose(), header=T)`

and copy and paste this into the R window at the command prompt; then press ENTER. This command opens a dialogue box and gives the user access to their filesystem. Select the file `calA.csv`. Then Click Open. The calibration data will then have been entered into the RStudio window.



6. Repeat step 5 to enter the buffer data (bufA.csv) and the input parameters (parA.csv) into RStudio.
7. Now the scene has been set and we can run the LOME program. Select the lines containing caout1, as shown below, and copy and paste them into the RStudio window, and press ENTER. The output will appear in the RStudio window. Also a plot of the fitted model and data will appear in the plot window, as shown in the next plot below.

LOMEscript.txt

```

1 # LOME for Fluorescence and Aequorin data, using R
2
3 # First, set your working directory by changing the next line. This tells
4 # R where your input files are.
5
6
7 setwd("~/Desktop/CaLOME")
8
9
10 # Now enter the program LOME
11
12 source("LOME.txt")
13
14 # Now select, one by one, three files that contain:
15
16 # the calibration data
17
18 # the buffer data
19
20 # the starting values for the parameters
21
22 cacal1 <- read.csv(file.choose(), header=T)
23
24 cabuf1 <- read.csv(file.choose(), header=T)
25
26 capar1 <- read.csv(file.choose(), header=T)
27
28 # Run the LOME program
29
30 caout1 <- lome(cacal1, cabuf1, capar1)
31
32 # Print the outputs. The outputs and a plot have been saved to
33 # your working directory that was set above.
34
35 caout1
36

```

RStudio


Console

```

~/Desktop/CaLOME/
> setwd("~/Desktop/CaLOME")
>
> source("LOME.txt")
> cacal1 <- read.csv(file.choose(), header=T)
> cabuf1 <- read.csv(file.choose(), header=T)
> capar1 <- read.csv(file.choose(), header=T)
> caout1 <- lome(cacal1, cabuf1, capar1)
>
> # Print the outputs. The outputs and a plot have been saved to
> # your working directory that was set above.
>
> caout1
[[1]]
  labels      parout
1 LigConc 9.000227e-01
2 Purity 9.000227e+01
3 pKapp 6.449712e+00
4 Kapp 3.550485e-07

[[2]]
      pX      X
1 5.044351 9.029200e-06
2 5.210704 6.155960e-06
3 5.381582 4.153538e-06
4 5.560015 2.754131e-06
5 5.756861 1.750408e-06
6 5.996519 1.008048e-06
7 6.353692 4.429020e-07
8 6.680679 2.086031e-07
9 6.994078 1.013730e-07
10 7.353031 4.435766e-08

```



8. The program will have written the output data and saved the plot to the working directory CaLOME in files with names starting with OUT, as shown below.

The screenshot shows the CaLOME directory with the following files and their sizes:

Name	Date Modified	Size
bufA.csv	Today at 19:53	234 bytes
bufF.csv	Today at 19:54	127 bytes
calA.csv	Today at 19:54	126 bytes
calF.csv	Today at 19:55	126 bytes
LOME.txt	Today at 20:18	3 KB
LOMEmcript.txt	Today at 22:05	741 bytes
OUTbufA.csv	Today at 22:03	388 bytes
OUTparA.csv	Today at 22:03	123 bytes
OUTplotA.pdf	Today at 22:03	8 KB
parA.csv	Today at 19:55	39 bytes
parF.csv	Today at 19:56	39 bytes